



ARTWORK  
PRINT DIRECTION

VARIABLE  
PRINT DIRECTION

Replace this box with GS1-128 barcode including production identifiers (if applicable) (01)50192253036631

LEGEND:		NOTES:		
	BREAK LINE		EYE MARK	ARTWORK PRINT PROCESS: THERMAL TRANSFER VARIABLE PRINT PROCESS: THERMAL TRANSFER  VARIABLE PRINT INFORMATION CONTAINS BARCODES, LOT NUMBER, EXPIRATION DATE, DATE OF MANUFACTURE AND SEQ. NUMBER WHERE APPLICABLE PER SKU.  MUST RESPECT SIZE, LOCATION AND ORIENTATION OF VARIABLE PRINT AREAS.  PRODUCTION IDENTIFIERS: SYMBOLS FOR LOT, USE-BY AND DoM ARE PART OF ARTWORK AND INCLUDED ONLY WHERE APPLICABLE PER SKU.
	CENTER LINE		ISO COLOR	
	CUT LINE		LABEL	
	DIMENSION		NO VARNISH	
	FOLD LINE		NO INK	
	1/8" INK BLEED		SEAL	
	PERFORATED LINE			
	PLANT NOTES			
	PRODUCT			
	TEAR LINE			
	TEXT			
	TITLE BLOCK			
	ZIPPER LINE			

**Instructions:**  
Barcodes for Cardinal Health products need to be generated in GS1-128, GS1 Data Matrix, UPC-A, or ITF-14 format. Barcodes must verify as an ANSI grade C (1.5) or better. The barcode must have human readable information below it (if there is a space constraint, the human readable information (HRI) may also be located around it, but clearly associated with the barcode). Human readable (HRI) dates must be in YYMMDD format. Additional (Non-HRI) date fields (e.g. "Use by:") should be ISO formatted as YYYY-MM-DD. For shipping cases, GS1-128 or ITF-14 are the only acceptable formats. For FDA Class I devices, the device identifier (DI) is the only information that is required to be encoded in the barcode. For FDA unclassified, FDA Class II, and FDA Class III devices, a minimum of one production identifier (PI) must be encoded in addition to the device identifier (DI). Although not a specification, additional guidance and information can be obtained in the Cardinal Health UDI Cookbook.

CardinalHealth™		Waukegan, IL 60085	
File Name:	8881201449_DL29A80144902_StandardNS_cs1bl	Colors:	Black
Dimensions:	8.000"W x 4.000"H		
Dieline Name:	5242101(8x4)-U_Rev2_PR		
Project #:	RB0020 Deland N&S No CE Mark Set 7		
Rev. Date:	Rev. 0   2021-05-14		
Designer:	Boss Group/Cella   C. Robison		
Comments/Printer Notes:	FPO barcode graphic represents the min acceptable size for accurate grade scan and cannot be reduced.		