

Scitara DLX™ Waters™ Empower™ 3 Connector

The Scitara DLX™ Waters™ Empower™ 3 Connector allows users to access the Empower application to create and update sample set methods and retrieve results. This Application Connector is installed directly onto a Windows-based computer.



Benefits

Enhanced Data Analysis

Utilizing data processing and analysis capabilities, the Waters Empower 3 Connector allows scientists to efficiently analyze chromatographic data, identify peaks, integrate signals, and perform quantification with high accuracy.

Improved Workflow

Creating Orchestrations using the Waters Empower 3 Connector streamlines the chromatography workflow by automating routine tasks such as data acquisition, processing, and reporting. This saves time and reduces the potential for human errors.

Customization and Flexibility

The Connector offers flexibility for customization according to specific laboratory requirements. Users can retrieve results in DLX based on specific sample sets or results, channel and injection IDs, and data processed or acquired.


Use Case

Waters™ Empower™ 3 in DLX

Create Empower Sample Set Method (SSM) & Retrieve Results

Executed by user

Executed by DLX



Working with CDS systems like Empower through the Scitara DLX Connector for Waters Empower 3 allows for bi-directional communication between LIMS or ELN and the Empower application. Users can send samples or sample sets to Empower, then retrieve and process results. The Waters Empower 3 Connector streamlines workflows, ensuring efficiency and accuracy.

For more information on how to build a connected laboratory:

[BOOK A DEMO HERE](#)