



Infusion Therapy Products



QUEST Medical, Inc.



Quest Medical, Inc. is committed to partnering with clinicians worldwide to produce better patient outcomes. We engineer infusion therapy products for superior infection control, organization, versatility and ease of use. Our innovative and high quality IV delivery sets are used to supply vital fluids and medications to patients in every department of your hospital or clinic.

Quest Medical aims to create efficiencies with each partnership through:
Inventory Management
Reliable Client Services
Next Day Shipping

Our representatives are subject matter experts who will assist you in standardizing products across departments, resulting in cost savings, reduced waste and faster set-up time.

To locate a sales specialist in your area, contact **800-672-0226** or email **custserv@questmedical.com**.

Icon Description Key



Contains Q2® Valve Technology



Low Volume (<3mL)



Contains Aspiration Port



European Conformity



Lead Time: 30 working days

Q2® NEEDLELESS CONNECTORS

3150
3150HF
3160

6

MANIFOLDS

95786
95787
95770
95789

7

95781
95782
95791
95788
95779
95778
95712

8

95711
9520
95710
95704
95702
95713
95707

9

HIGH-FLOW STOPCOCKS

HFSTV01
HFSTV02
HFSTV03
STV01
STV02
STV03
ST01

10

ST02
ST03

11

PRECISION DELIVERY® EXTENSION SETS

PD06F
PDM06F
PDMV06F
PDV06
PDV06F

12

PD60
PD60F
PDM60
PDM60F
PDMV60
PDMV60F
PDV60

13

PDV60F
PDBV05
PDB05F
PDB05DF
PDBV05F
PDBV05DF

14

PDT05
PDT05F
PDT05DF
PDTV05
PDTV05F
PDTV05DF

15

LIGHT-PROTECTIVE, AMBER EXTENSION SETS

PDV06A
PDV60A
PDMV60FA
PDV60FA
PDTV05A
PDTV05FA

16

PDTV05DFA

17

EXTENSION SETS

95701
9219
95902
95725S
95724

18

95725
95726
95727
LOW-PRESSURE, POWER
INJECTION EXTENSION SETS
95714
95906
95907

19

IV ADMINISTRATION SETS

22-201-V
ASV101Y
ASV102Y
ASV103Y

20

ASV60125
ASV10125
ASV10115
ASV10130
95780

21

MICROBORE ADMINISTRATION SETS

ASV170F
ASV370
ASV470

22

Our Evolving Product Lines


Quest Medical, Inc. was founded on the premise of developing medical devices that would have a positive impact on patient care. Since 1979, we have been dedicated to providing the highest quality products to our three focused markets: Infusion Therapy, Cardiovascular and Ophthalmology.

Quest's evolving Infusion Therapy Line features rigid and flexible manifolds with removable components and intravenous (IV) administration sets with high or low volume configurations to fit all patients' needs. Our products incorporate the following trademarks: Q2® Valve Technology (Needleless Connectors/NCs), Precision Delivery® Extension Sets, ExoFlex® Manifolds, MultiPort® Manifolds and Qlinch™ Securement Bands.



Integrated into most of Quest's Infusion Products is our Q2® Valve Technology. Gaining access to each patient's bloodstream via Infusion Therapy can increase risk of infection. Minimize that risk by incorporating a simple, closed-system design into your line setup. Quest's Q2® Valve Technology is void of internal componentry, reducing the chance of contamination.

“When needleless systems are used, a split septum valve may be preferred over some mechanical valves due to increased risk of infection with the mechanical valves.”



- CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services



Understanding patient safety is your #1 priority, Quest has revolutionized needleless connectors with our Q2[®] Valve Technology. Integrated into most Quest products, the enhanced safety features of Q2[®] NCs mitigate infection risks. The straight, unobstructed fluid path of Q2[®] NCs provides high flow rates and is engineered to suppress thrombotic occlusion risks caused by disconnecting syringes and other sets. Q2[®] NCs enable multiple activations and maintain a sterile, closed system for increased patient safety and improved outcomes, providing you with peace of mind.



3150	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2 [®] Needleless Connector, swabable	N/A	N/A	0.09 mL	200



3150HF	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2 [®] High-Flow Needleless Connector, swabable		N/A	N/A	0.09 mL	200




3160	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2 [®] Vial Access Adapter, swabable	N/A	N/A	0.06 mL	50




To help clinicians accommodate an ever-changing patient population, Quest offers safe and simple solutions for fluid management. From high-flow stopcocks to the rigid MultiPort[®] Manifold and flexible ExoFlex[™] Manifold, Quest designed the complete manifold line to adapt to your needs. The closed-system, procedure-ready sets aid in decreasing line manipulation and lowering infection risks, providing you even more time to devote to patient care.




95786	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex [™] Manifold with Q2 [®] NC + integrated check valve, (4) double T-sites with Q2 [®] NCs + integrated check valves, removable spines and slide clamp	0.06 in.	11.25 in. (29cm)	1.75 mL	50




95787	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex [™] Manifold with Q2 [®] NC + integrated check valve, (6) T-sites with Q2 [®] NCs + integrated check valves, removable spines and pinch clamp	0.06 in.	14.75 in. (37.5cm)	1.65 mL	50




95770	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex [™] Manifold with Q2 [®] NC + integrated check valve, (2) double T-sites with Q2 [®] NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	9 in. (22cm)	1.2 mL	50



95789	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex [™] Manifold with Q2 [®] NC + integrated check valve, (2) double T-sites with Q2 [®] NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2 [®] NC	0.10 in.	53 in. (135cm)	6.6 mL	50




95781	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	10.5 in. (26.75cm)	1.75 mL	50
				


95782	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	10.5 in. (26.75cm)	1.17 mL	50
				


95791	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	16 in. (40.65cm)	3 mL	50
				


95788	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	55 in. (140cm)	7.3 mL	50
				


95779	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	53 in. (135cm)	6.7 mL	50
				


95778	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	9 in. (22cm)	1.2 mL	50
				


95712	TUBING ID	APPROX. LENGTH	APV	CASE QTY
MultiPort® Manifold with (6) Q2® NCs + integrated check valves, double T-site with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	2.2 mL	50
				


95711	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (3) double T-sites with Q2® NCs + integrated check valves, (2) slide clamps and additional T-site with Q2® NC + integrated check valve	0.06 in.	8.5 in. (22cm)	1.5 mL	50
				

9520	TUBING ID	APPROX. LENGTH	APV	CASE QTY
MultiPort® Manifold with (6) Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	2.0 mL	50
				

95710	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (6) T-sites with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	14.5 in. (37cm)	1.0 mL	50
				

95704	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves and slide clamp	0.06 in.	7 in. (18cm)	1.1 mL	50
				

95702	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	0.9 mL	50
				

95713	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves and side clamp (Not shown: T-sites enclosed in green casing)	0.06 in.	10 in. (25cm)	0.9 mL	50
				

95707	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs + integrated check valves	0.06 in.	4 in. (10cm)	0.7 mL	50
				

High-Flow Stopcocks

HFSTV01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® High-Flow 4-Way Stopcock Extension Set with Q2® NC, with pinch clamp + Y-site with Q2® NC	0.12 in.	43.8 in. (111cm)	7.7 mL	50





HFSTV02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon Q2® 2-Gang High-Flow 4-Way Stopcock Infusion Set with Q2® NC, with pinch clamp + Y-site with Q2® NC	0.12 in.	45.5 in. (115cm)	8.0 mL	50





HFSTV03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon Q2® 3-Gang High-Flow 4-Way Stopcock Infusion Set with Q2® NC, with pinch clamp + Y-site with Q2® NC	0.12 in.	47.1 in. (119cm)	8.3 mL	50





STV01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	.5 mL	100





STV02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon Q2® 2-Gang High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	0.9 mL	50





STV03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon Q2® 3-Gang High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	1.1 mL	50





ST01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon High-flow 4-Way Stopcock	N/A	N/A	0.6 mL	100





High-Flow Stopcocks continued

ST02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon 2-Gang High-Flow 4-Way Stopcock	N/A	N/A	0.8 mL	50





ST03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon 3-Gang High-Flow 4-Way Stopcock	N/A	N/A	1.5 mL	50









Quest Medical's Precision Delivery® Infusion Therapy Line is intentionally engineered to provide safe and ultra-precise infusions for your youngest critical-care patients. Specialized componentry minimizes dilution and priming volume while delivering all solutions to neonatal and pediatric patients. Reduce the risk of CLABSI's using our Q2® Valve Technology. A variety of configurations with specialized componentry such as ultra-microbore tubing, color-coded slide clamps, light-protective, amber tubing and pediatric-specific filters, provide for ultimate control of infusions.




PD06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 1.2 micron filter	0.04 in.	7.8 in. (19cm)	0.4 mL	50



PDM06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 0.2 micron filter	0.02 in.	7.8 in. (19cm)	0.3 mL	50




PDMV06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and 0.2 micron filter	0.02 in.	8 in. (20cm)	0.3 mL	50




PDV06	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC and slide clamp	0.04 in.	7.75 in. (19cm)	0.4 mL	50




PDV06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and 1.2 micron filter	0.04 in.	8.5 in. (21cm)	0.7 mL	50



PD60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp	0.04 in.	61 in. (155cm)	1.4 mL	100



PD60F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 1.2 micron filter	0.04 in.	62 in. (158cm)	1.7 mL	100



PDM60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and ultra-microbore tubing	0.02 in.	61 in. (155cm)	0.5 mL	100




PDM60F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp, 0.2 micron filter and ultra-microbore tubing	0.02 in.	62 in. (158cm)	0.6 mL	100



PDMV60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and ultra-microbore tubing	0.02 in.	61 in. (155cm)	0.6 mL	100




PDMV60F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp, 0.2 micron filter and ultra-microbore tubing	0.02 in.	62 in. (158cm)	0.8 mL	100




PDV60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC and slide clamp	0.04 in.	61 in. (155cm)	1.6 mL	100




PDV60F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and 1.2 micron filter	0.04 in.	62 in. (158cm)	1.6 mL	100




PDBV05	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps and 2-in-1 male luer	0.04 in.	6 in. (16.5cm)	0.33 mL (W) 0.33 mL (B)	50




PDB05F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Bifuse Extension Set with (2) slide clamps, 0.2 micron filter and 2-in-1 male luer	0.04 in.	7.28 in. (18.5cm)	0.38 mL (W) 0.55 mL (Y)	50




PDB05DF	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Bifuse Extension Set with (2) slide clamps, 0.2 micron filter, 1.2 micron filter and 2-in-1 male luer	0.04 in.	7.28 in. (18.5cm)	0.55 mL (W) 0.55 mL (Y)	50




PDBV05F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps, 0.2 micron filter and 2-in-1 male luer	0.04 in.	7.4 in. (19cm)	0.33 mL (W) 0.50 mL (Y)	50




PDBV05DF	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps, 0.2 micron filter, 1.2 micron filter and 2-in-1 male luer	0.04 in.	7.4 in. (19cm)	0.50 mL (W) 0.50 mL (Y)	50




PDT05	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Trifuse Extension Set with (3) slide clamps and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.2 mL (W) 0.2 mL (B) 0.2 mL (Y)	50




PDT05F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Trifuse Extension Set with (3) slide clamps, 0.2 micron filter and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.2 mL (W) 0.2 mL (B) 0.4 mL (Y)	50




PDT05DF	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Trifuse Extension Set with (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.4 mL (W) 0.2 mL (B) 0.4 mL (Y)	50




PDTV05	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.25 mL (Y)	50



PDTV05F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps, 0.2 micron filter and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.50 mL (Y)	50



PDTV05DF	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.50 mL (W) 0.25 mL (B) 0.50 mL (Y)	50



Light-Protective, Amber Extension Sets

PDV06A	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC and slide clamp	0.04 in.	7.75 in. (20cm)	0.4 mL	50




PDV60A	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC and silde clamp	0.04 in.	61 in. (155cm)	1.5 mL	50




PDMV60FA	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with ultra-microbore, light-protective, amber tubing, Q2® NC, slide clamp and 0.2 micron filter		0.02 in.	62 in. (158cm)	0.8 mL	50




PDV60FA	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC, slide clamp and 1.2 micron filter	0.04 in.	62 in. (158cm)	1.6 mL	50



PDTV05A	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.25 mL (Y)	50




PDTV05FA	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps, 0.2 micron filter and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.50 mL (Y)	50



Light-Protective, Amber Extension Sets continued


PDTV05DFA	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps, 0.2 micron filter and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.50 mL (W) 0.25 mL (B) 0.50 mL (Y)	50



Quest's Extension Sets focus on safe and accurate infusion therapy delivery with minimal priming volumes and high flow rates. Our offerings include various tubing sizes and component options including low-pressure, power injectable extensions that allow you to connect directly to your patient's catheter when administering contrast for diagnostic imaging.



95701	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Extension set with Q2® NC + integrated check valve, (1) T-site with Q2® NC + integrated check valve and pinch clamp	0.06 in.	7 in. (18cm)	0.66 mL	50




9219	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Female Y-extension set with (1) female luer pressure-activated check valve, (1) Q2® NC + integrated check valve and pinch clamp	0.10 in.	42 in. (107cm)	6.2 mL	50



95902	TUBING ID	APPROX. LENGTH	APV	CASE QTY
High-flow filter extension set with high-flow 0.22 micron filter, pinch clamp and Y-site with Q2® NC	0.10 in.	17 in. (43.2cm)	4.5 mL	100



95725S	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with 24 in. perforated sleeve and slide clamp	0.045 in.	25 ft. (7.6m)	10 mL	50



95724	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Extension set with pinch clamp	0.10 in.	20 ft. (6m)	31 mL	100



95725	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	20 ft. (6m)	8.5 mL	100



95726	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	10 ft. (3m)	3.2 mL	50




95727	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	60 in. (152cm)	1.7 mL	100




Low-Pressure, Power Injection Extension Sets
Up to 325 psi


95714	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with detachable Q2®NC and pinch clamp	0.10 in.	7 in. (18cm)	1.0 mL	100



95906	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with slide clamp and Q2® NC + T-connector extension	0.045 in.	6 in. (15cm)	0.3 mL	100




95907	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with Q2® NC side port, slide clamp and Q2® NC + T-connector extension	0.045 in.	6 in. (15cm)	0.4 mL	100




Quest's IV Administration Sets assist with accurate delivery of medication. Accessing the patient's bloodstream when administering drugs may increase the risk of infection. With the incorporation of Q2® Needleless Connectors, Quest's IV Administration Sets are engineered with a closed system to promote patient safety. Multiple tubing ID options and detachable extensions provide a variety of configurations, which meet the demands of various departments in your facility including Anesthesia, Oncology and Critical Care.




22-201-V	TUBING ID	APPROX. LENGTH	APV	CASE QTY
20-drop/mL vented secondary set with ball shut-off valve and roller clamp	0.10 in.	75 in. (190cm)	15.5 mL	48




ASV101Y	TUBING ID	APPROX. LENGTH	APV	CASE QTY
10-drop/mL administration set with check valve, (1) Y-site connectorn + Q2® NC, roller clamp, (1) T-site Q2® NC + integrated check valve and pinch clamp	0.10 in.	90 in. (227cm)	15 mL	50




ASV102Y	TUBING ID	APPROX. LENGTH	APV	CASE QTY
10-drop/mL administration set with check valve, (2) Y-sites + Q2® NCs, roller clamp, (1) T-site with Q2® NC and pinch clamp	0.10 in.	99 in. (251cm)	18 mL	50




ASV103Y	TUBING ID	APPROX. LENGTH	APV	CASE QTY
10-drop/mL administration set with check valve, (3) Y-site with Q2® NCs, roller clamp, (2) T-sites with Q2® NCs + integrated check valve and pinch clamp	0.10 in.	107 in. (272cm)	20 mL	50




ASV60125	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
60-drop/mL administration set with check valve, (4) Y-sites with Q2® NCs, (2) pinch clamps, roller clamp, detachable high-flow 4-way stopcock + Q2® NC and detachable 6.5 in. extension with Q2® NC		0.10 in.	126 in. (320cm)	29.1 mL	50




ASV10125	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, detachable high-flow 4-way 2-gang stopcock + Q2® NCs and pinch clamp		0.10 in.	129 in. (327cm)	28.6 mL	50




ASV10115	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
10-drop/mL administration set with check valve, Y-site with Q2® NC, roller clamp, detachable high-flow 4-way stopcock + Q2® NCs and detachable high-flow 4-way 2-gang stopcock + Q2® NCs		0.10 in.	115 in. (292cm)	27.9 mL	50



ASV10130	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, detachable high-flow 4-way 3-gang stopcock + Q2® NCs, with detachable 6 in. extension with Q2® NC and pinch clamp		0.10 in.	132 in. (335cm)	30.4 mL	50




95780		TUBING ID	APPROX. LENGTH	APV	CASE QTY
10-drop/mL administration set with check valve, Y-site + Q2® NC, roller clamp, (3) T-sites with Q2® NCs + integrated check valves with ExoFlex™ Manifold + removable spine and detachable extension set with (2) Y-sites + Q2® NCs and pinch clamp		0.10 in.	117 in. (297.18cm)	22 mL	50




Microbore Administration Sets


ASV170F	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore administration set with high-flow spike, pressure activated valve, pinch clamp, Q2® NC syringe port, check valve, slide clamp and 0.2 micron filter		0.10/0.04 in.	107 in. (271cm)	8.5 mL	50



ASV370	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore trifurcated administration set with high-flow spike, pressure activated valve, pinch clamp, (3) Q2® NC syringe ports, check valve and (3) slide clamps		0.10/0.04 in.	103 in. (261cm)	0.3 mL (W) 0.3 mL (Y) 8.9 mL (B)	50



ASV470	Coming Soon	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore quadfurcated administration set with high-flow spike, pressure activated valve, pinch clamp, (4) Q2® NC syringe ports, check valve and (4) slide clamps		0.10/0.04 in.	103 in. (261cm)	0.3 mL (R) 0.3 mL (W) 0.3 mL (Y) 8.9 mL (B)	50



Ordering Process and Terms

To place an order for Infusion Therapy Products, please see below for ordering process and terms.

Distribution/Order Placement

Quest Medical, Inc.
One Allentown Parkway
Allen, Texas 75002

Phone: 800.627.0226 x 461
Fax: 972.390.7173
Email: Custserv@questmedical.com

Vendor Federal Tax I.D. No.
75-274056

F.O.B. ORIGIN
F.O.B., Pre-pay, and Add

Minimum Order
Need to order in case quantity

Delivery Lead Time
Lead time will be no more than forty-eight (48) to seventy-two (72) hours following receipt of order.
*Unless specified with 30-day lead time

Payment Terms
Net 30 days, with approval

Remit Address
Quest Medical, Inc.
PO Box 934486
Atlanta, GA 31193-4486

A/R Phone: 800.627.0226 x 271 or 299

*Custom Setup Inquiries contact
Samples@questmedical.com

CONNECT WITH QUEST!





QUEST Medical, Inc.

custserv@questmedical.com | +1 800.627.0226
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com



EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Rx ONLY  Made in USA