



BKACTIV

bk medical 
a GE Healthcare company

NEXT LEVEL DECISION SUPPORT

Real-time intraoperative ultrasound for neurosurgery and spine

EMPOWERING SMARTER SURGERY

bkActiv is designed for all neurosurgical subspecialties and offers features, such as color Doppler and compare modes, that enhance surgical workflow in even more procedures.

Our pioneering visualization and guidance platform delivers high-resolution images and data-driven insights when it matters most — during the procedure.




UNLOCK REAL-TIME IMAGING FOR MORE PROCEDURES

NEURO-ONCOLOGY

MINIMALLY INVASIVE



Metastatic lesion,
Dual Live Compare.
N13C5 Transducer



Pituitary Adenoma adjacent
to Carotid Artery.
N20P6* Transducer



NEUROVASCULAR



Hematoma due to AVM
with color Doppler.
N13C5 Transducer



SPINE



Spinal Decompression,
Dual Live Compare.
X18L5s** Transducer



PEDIATRICS



Shunt Placement into Brain
Cyst, Pediatrics.
N11C5s *** Transducer



DESIGNED FOR THE NEUROSURGEON

bkActiv is real-time active imaging that takes intraoperative ultrasound to the next level and is tailored to the neurosurgical operating room.



Support image interpretation with review and compare modes



Extended reach monitor arm fits into small, crowded ORs. Switch between the monitor's portrait and landscape options.

Customize to your needs with a touchscreen interface that can be adapted to suit your OR staff.



Use the sterilizable remote for even more control.



TRANSDUCERS



N20P6*
Transducer for hard-to-reach areas with small openings.



X18L5s**
Transducer for extreme high resolution.



N13C5
Transducer with wide field of view for full visibility.



N11C5s***
Transducer with deep penetration and small footprint.

REAL-TIME VISUALIZATION HAS NEVER LOOKED THIS GOOD

TAKING IMAGE QUALITY EVEN FURTHER

Prism Technology powers bkActiv with unique algorithms like Tru-Focus that are specifically designed for imaging the brain and spine.



COMPARE MODES

Support image interpretation during your procedures with three different compare features that enable you to compare live imaging with previously taken scans.



CUSTOMIZABLE USER INTERFACE

A simplified user interface to see only what you need. Customize by adding or removing controls.



*The N20P6 Transducer has not been licensed by Health Canada or CE-marked.
**Use of the X18L5s Transducer for intraoperative (neuro) has not been CE-marked.
***The N11C5s Transducer has been cleared for adults and children (>2 - 12 years).

USA
Sales & Service
25 Corporate Drive,
Suite 230, Burlington,
MA 01803, USA
T +1 978 326 1300
F +1 978 326 1399

EUROPE AND REST OF WORLD
Sales, Service & Design Center
Mileparken 34,
2730 Herlev,
Denmark
T +45 4452 8100
F +45 4452 8199

BK MEDICAL APS
Manufacturer
Mileparken 34,
2730 Herlev,
Denmark
T +45 4452 8100
F +45 4452 8199